## Calibrating Unit New Standard Solution Report

Equipment Location: Calibration File No.: Certification File No.: Linearity File No.: Solution File No.: Sequential File No.:	Alcotest 7110 I CLIFTON CIT 01040 01041 01042 01052 01052		Calib. Date: Cert. Date: Lin. Date: Soln. Date: File Date:	02/15/2010 02/15/2010 02/15/2010 02/15/2010 03/08/2010 03/08/2010	Serial No.: Calib. No.: Cert. No.: Lin. No.: Soln. No.:	00011 00011
Calibrating Unit: Control Solution %: Solution Control Lot:	WET 0.100% 09A062		Model No.:	CU-34	Serial No.: Expires: Bottle No.:	DDXA \$3-0032 01/15/2011 1184
Function		Result	Time	Temperature	Com	ment(s)
Tunction		%BAC	HH:MM	Simulator (°C)	or E	rror(s)
Ambient Air Blank		0.000%	10:09S			- Lawrence tolet
Control 1 EC		0.102%	10:10S	34.0°C		PASSED ***
Control 1 IR		0.101%	10:10S	34.0°C	*** TEST	PASSED ***
Ambient Air Blank		0.000%	10:10S			
Control 2 EC		0.100%	10:11S	34.0°C		PASSED ***
Control 2 IR		0.101%	10:11S	34.0°C	*** TEST	PASSED ***
Ambient Air Blank		0.000%	10:12S			
Control 3 EC		0.101%	10:12S	34.0°C		PASSED ***
Control 3 IR		0.101%	10:12S	34.0°C	*** TEST	PASSED ***
Ambient Air Blank		0.000%	10:13S			*

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in acordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: BRA

Signatura:

First Name QAVID

MI: A

Badge No.: 10

Date:

03/08/2010



JON S. CORZINE Governor

## State of New Jersey

OFFICE OF THE ATTORNEY GENERAL DEPARTMENT OF LAW AND PUBLIC SAFETY DIVISION OF STATE POLICE POST OFFICE Box 7068 West Trenton NJ 08628-0068 (609) 882-2000

ANNE MILGRAM Attorney General

COLONEL JOSEPH R. FUENTES Superintendent

## CERTIFICATION OF ANALYSIS . 0.10 PERCENT BREATH ALCOHOL SIMULATOR SOLUTION

ACCEPTANCE SPECIFICATIONS FOR BREATH ALCOHOL SIMULATOR SOLUTION: Ethyl alcohol concentration within, but not exceeding, the range of 0.117 to 0.125 grams per 100 milliliters of solution.

MANUFACTURER: Drager Safety, Inc.

**ANALYSIS DATE: 2/11/2009** 

## BREATH ALCOHOL SIMULATOR SOLUTION LOT NUMBER: 09A062

Representative samples of the above-referenced Lot Number were tested by Gas Chromatography and found to have an ethyl alcohol concentration range of 0.1195 to 0.1197 grams per 100 milliliters of solution.

This lot of breath alcohol simulator solution may be utilized as a known traceable standard for the purpose of conducting periodic tests, pursuant to N.J.A.C. 13:51-3.4, of approved breath test instruments (N.J.A.C. 13:51-3.5) utilized by law enforcement agencies in this State. The manufacturer's expiration date for this lot of breath alcohol simulator solution is January 15, 2011.

As Chief Forensic Scientist of the Division of State Police, I hereby certify and attest that the tests and results documented in this Certificate of Analysis were performed at my direction and under my supervision by personnel of, and at, the Office of Forensic Sciences of the Division of State Police on properly functioning and calibrated instruments and equipment. All procedures utilized are accurate, objective, and performed on a routine basis by personnel within the Office of Forensic Sciences, in accordance with their professional duties and responsibilities.


Chief Forensic Scientist

Division of State Police

Sworn to and subscribed before me this of O day of Juliu

Linda L. DeSantis My Commission

Expires Aug. 17, 2009



DEPARTMENT OF

THE AND FILLIE SEFET

ODIVIDE BEARCAGE

CHIEF TO COMMETCH TO CO

RIGINAL COU	RSE DATES	3-29-2006			
DATE	Refresher Course	INSTRUCTOR			
2-7-08	Sussex Co. FA	Russo			
	,				
•					
i.					
i,					
7.		•			
t.					
9.					